

U95

1.95Litre Ultrasonic Bath

Ultrawave's **U-Series** is compact and user friendly, offering effective cleaning and processing of small items in both professional and domestic applications.

These include healthcare, podiatry, tattooists, jewellers and many more.

Available in heated and unheated versions, the **U-Series** allows simple control of time (0 to 15 minutes) and temperature (ambient to 70°C) where heated options are selected.



The **U95** is a 1.95Litre ultrasonic bath with a timer control of up to 15 minutes.

Each **U-Series** is supplied complete with a stainless steel basket and ABS plastic lid, useable as a drip tray to collect excess water when the basket is removed from the tank.

The **U95** provides the cost benefits of excellent entry level ultrasonics with Ultrawave's long renowned quality and reliability.



U95

1.95Litre Ultrasonic Bath

Construction Case	304 Stainless Steel
Construction Tank	304 Stainless Steel
Control	Analogue (Time)
Overall Dimensions	210 x 165 x 160mm (WxDxH)
Tank Dimensions	200 x 150 x 75mm (WxDxH)
Basket Dimensions	180 x 135 x 45mm (WxDxH) (60mm height with handles)
Maximum Capacity	2.25Litres
Working Capacity	1.95Litres
Ultrasonic Power	35W
Ultrasonic Power (per litre)	17.9W/litre
Operating Frequency	44KHz
Maximum Heating Capacity	N/A
Heater Power	N/A
Timer	0 - 15 minutes
Mains Connection	220 – 240V (50 - 60Hz)
Drain Outlet	N/A
Lid	ABS Plastic
Basket	Stainless Steel
Gross Weight (inc. lid, basket & packaging)	4.5KGs
Net Weight	3.7KGs
Warranty	3 Years

